

Features	Z-ONE	ZC3	ZH3	ZT-Clinic MED	NICE
plug & play - no preparation needed	Yes	Yes	Yes	Yes	prepare water before use
Destination of use	COLD	Gradual Contrast - Multi Phase Temperature (H/C)	Gradual Contrast - Multi Phase Temperature (H/C)	Gradual Contrast - Multi Phase Temperature (H/C)	COLD
Thermal System	Heat Pump Refrigeration System	Heat Pump Refrigeration System	Heat Pump Refrigeration System	Heat Pump Refrigeration System	Peltier System
Heat Extraction Power	450w	750w	750w	1500w	150w
Heat exchange fluid	Non Tox Freeze (Glycol)	Non Tox Freeze (Glycol)	Non Tox Freeze (Glycol)	Non Tox Freeze (Glycol)	Water
Range temperature	All Temp° Between 5-15°C	All Temp° Between 0-40°C	All Temp° Between 0-40°C	All Temp° Between -5 - 40°C	5 Levels from 5,5 to 14,4°C
Preloaded System Liquid	Yes	Yes	Yes	Yes	No
liquid connectors	Dripless system	Dripless system	Dripless system	Dripless system	Leaky system
air compression	Yes	Yes	Yes	Yes	Yes
closed circuit system	Yes	Yes	Yes	Yes	No
Contemporary Wraps	up to 2	up to 2	up to 2	up to 4	Only 1
preloaded Hoses	Yes	Yes	Yes	Yes	No
preloaded wraps	Yes	Yes	Yes	Yes	No
anatomic Wraps	Yes	Yes	Yes	Yes	Yes
Air compression Levels	40-60-80 mmHg	40-60-80 mmHg	40-60-80 mmHg	40-60-80mmHg	13-26-39 mmHg
Active Air Compression Duration settings	3" / 1' / 2' / 3' / constant	3" / 1' / 2' / 3' / constant	3" / 1' / 2' / 3' / constant	3" / 1' / 2' / 3' / constant	only 2' ON
weight	9kg	14kg	26kg	48kg	5kg (empty)
dimension	cm 25 x 28 x 30	cm. 32 x 38 x 33h	cm. 21 x 26 x 86	cm 36 x 55 x 93	cm 22 x 22 x 28
Maintenance	No maintenance	No maintenance	No maintenance	No maintenance	Drain off all the liquid inside the wraps and the unit every time you complete a treatment